

AIR-FLOW®

3.0




AIR-FLOW® handy 3.0

EMS⁺
ELECTRO MEDICAL SYSTEMS

AIR-FLOW® handy 3.0
ORIGINAL AIR-FLOW® THERAPY

AIR-FLOW® handy 3.0

ORIGINAL AIR-FLOW® THERAPY WITH A NEW STANDARD IN EFFICIENCY
AND FLEXIBILITY, ERGONOMICS AND DESIGN



PREMIUM PROPHYLAXIS 3.0

MOBILE PROPHYLAXIS DEVICE MEETING THE MOST
STRINGENT REQUIREMENTS OF MODERN
PROFESSIONAL TOOTH CLEANING OVER THE ENTIRE
TREATMENT SPECTRUM BOTH **SUB- AND SUPRAGINGIVAL**


MODERN PROPHYLAXIS

NEW OPTIONS FOR SUCCESSFUL THERAPY →



AIR-FLOW® handy 3.0 PREMIUM

> AIR-FLOW handy 3.0 PLUS with additional AIR-FLOW PERIO handpiece – always ready for sub- and supragingival use with the AIR-FLOW® Powders PERIO and PLUS




AIR-FLOW® handy 3.0 PLUS

> Optimal for use in periodontal pockets up to 4 mm in depth as well as over the entire visible tooth area – flexible and efficient with the AIR-FLOW® Powder PLUS



AIR-FLOW® handy 3.0 PERIO

> For optimal subgingival use in pockets deeper than 4 mm – simple and thorough removal of biofilm with the AIR-FLOW® Powders PERIO and PLUS



AIR-FLOW® handy 3.0

> Specialized in all supragingival applications – for optimal results on the entire visible tooth area with the AIR-FLOW® Powders CLASSIC, CLASSIC COMFORT and SOFT



BEAUTIFUL → BEAUTIFULLY LIGHT

DESIGN AND FUNCTIONALITY FOR OPTIMUM TREATMENT COMFORT




> Designed to ensure easy vision of the treatment area

> Well-balanced for comfortable handling



AIR-FLOW® handy 3.0



**ORIGINAL AIR-FLOW®
PERIO THERAPY
PENETRATES INTO THE
DEEPEST PERIODONTAL
POCKETS – ELIMINATES
BIOFILM AND WITH IT
MILLIONS OF POTENTIAL
PATHOGENIC GERMS**

PERIO nozzle for subgingival application – three horizontal nozzle outlets for powder-air mixture, one vertical nozzle outlet for water

DENTAL IMPLANTS → PROPHYLAXIS WITH THE ORIGINAL AIR-FLOW® THERAPY AND THE AIR-FLOW® handy 3.0 PERIO




> Implant exposed due to bone loss

- > Peri-implantitis is the periodontitis of implant patients
- > Plaque adheres to the implant and causes inflammation around the implant, which in turn causes bone loss and loosening of the implant
- > Regular subgingival prophylaxis using the Original AIR-FLOW® Therapy prevents peri-implantitis



**AIR-FLOW®
EASY-FILL**
> Safe,
quick and
clean



FULL-MOUTH AIR-FLOW® THERAPY

A MIX OF AIR AND ORIGINAL AIR-FLOW®
POWDER* COMBINES WITH A JET OF WATER FOR
RAPID AND PAINLESS REMOVAL OF DENTAL PLAQUE
AS WELL AS BIOFILM FROM THE ENTIRE MOUTH –
MINIMALLY ABRASIVE AND DUST FREE

*All Original AIR-FLOW® Prophylaxis Powders are developed in EMS internal laboratories and are produced by EMS – for optimal performance with Original AIR-FLOW® handpieces



EMS-SWISSQUALITY.COM

EMS
ELECTRO MEDICAL SYSTEMS SA
Chemin de la Vuarpillière 31
CH-1260 Nyon

Tel. +41 22 99 44 700
Fax +41 22 99 44 701
welcome@ems-ch.com
www.ems-dent.com