



Features:

- Durable Zinc-alloy shell
- Advanced DSLR technology
- Special adjustable focus
- Resolution as 700 TV lines
- Less than 5% distortion

ProMi

Adjustable Focus Intraoral Camera HDI-712B


Working Scheme:



1.High-definition image quality with the distortion less than 5% could display the odontoschim perfectly.




2.Zero distance photography makes it easier to obtain fantastic images of posterior teeth for the patients with restriction problem of mouth opening .



3.Imaging efficacy close to those from dental operating microscope, which could observe the results of root canal cleaning and the root canal orifice after odontotrypy.



4.Adjustable focus which works like DSLR could fulfill the different image requirements from root canal to whole face.



ProMi



Specification

Sensor	1/3"HD CMOS
Effective Pixels	1280(H) × 720 (V)@30f
Resolution	700 TV lines
Angle of View	63°
Focus	1mm~infinity
Lighting	8 leds
Distortion	<5%
Pixel	1.0Mega pixels
Operation System	Windows2000/xp/vista/Win7/Win 8(32bit&64bit)

HANDY(UK/SHANGHAI)TECHNOLOGIES CO., LTD

Tel: +86 021-56766108 Fax: +86 021-56766308

Website: www.handydental.cn www.handycreate.com

E-mail: h.sale@hotmail.com handy@handycreate.com

After - sales supports:

Skype: [handysupports](https://www.skype.com/handysupports) Email: support@handycreate.com