



Features:

- Modern aluminum shell design
- Revolutionary iTouch button
- Perfect 720P HD images
- Super depth of field lens

TouchMi

Auto-focus (fixed) Intraoral Camera HDI-220B

Working Scheme:



Features:

- Revolutionary **iTouch** button brings brand new customer experience
- Modern **aluminum** shell design makes it concise and durable
- **Lens** with super depth of field make the view much clearer
- High resolution as **720P** for HD images



Will you still use the interoral camera with 640x480 resolution while the monitors over 20" are under popularization?

Specification	
Sensor	1/4 " 2 Mega Pixel CMOS
Resolution	640X480~1280X720@720P
Frame Rate	30fps@640x480 ~ 1280x 720@720P
Signal Output	USB 2.0
View of Angle	90°
Lighting	6 LED with High Luminance
Length of Cable	2.5m
Focus Range	5mm ~ 35mm
Operation System	Window 2000/XP/ Vista /Win7/Win8 (32bit & 64bit)

HANDY(UK/SHANGHAI)TECHNOLOGIES CO., LTD

Tel: +86 021-56766108 Fax: +86 021-56766308
Website: www.handydental.cn www.handycreate.com
E-mail: h.sale@hotmail.com handy@handycreate.com
After - sales supports:
Skytel: handysupports Email: support@handycreate.com