# Superiority is at hand – all day, every day. Maximum Power. Whisper Quiet. Optimal Comfort.

NEW KaVo MASTERtorque<sup>™</sup> with Direct Stop Technology

KaVo Master Series

A true master knows no compromise.

MASTER torque

X

MASTER forque

R torque M

MASTE

Silence Technology Small diameter Small head height No retraction



### KaVo. Dental Excellence.

## Maximum Power.

KaVo's MASTERtorque – experience a completely new dimension of power with 23 watts peakpower. Best conditions for comfortable and efficient preparation work.

### More power 23 watts peak power for approx. 20 % more torque\* Higher precision


Stable operating speed and high torque all the way through to the highest power level

More comfort Efficient and effective preparation work even at low contact pressure

\* Compared to KaVo's GENTLEsilence 8000 B

# Whisper Quiet.

KaVo's new MASTERtorque combines the proven Silence Technology with further frequency optimisation. Best conditions for quiet and comfortable preparation work even at full speed and during stopping of the turbine.



and low noise levels of only 57 dB(A)





2

2 torque

0

Ľ

Ш H ហ

∢

Σ

2

3

Maximum power:

for stable operating speed

# Relaxed work with Direct Stop Technology.


Rapid stop – for safety. The new Direct Stop Technology (DST) reduces the stopping time of KaVo's MASTERtorque by approx. 60 % and

\*\* on average

### Safe preparations

Innovative Direct Stop Technology stops the bur in one second\*\* and thus ensures increased safety during treatments

### Retractive suction from the head



stops in only one second\*\*. Best conditions for completely safe and relaxed work.








KaVo. Dental Excellence.

# Excellent visibility.

KaVo MASTERtorque's compact and ergonomic design together with KaVo's unique lighting concept create best working and visibility conditions at the preparation area.



### Better visibility

- Reduced head diameter on the bur side. Head height of only 33/64 in (13.1mm); with 3/4 in (19mm) dental bur only 51/64 in (20mm).
- With patented KaVo 100° / 19° head / knee angle combination

#### Head height in comparison to other manufacturers

- More light
  - Preparation areas optimally illuminated at all times
  - Preparation areas optimally visible in spray mode, light and spray nozzles are offset



Optimum visibility through offset spray and light exit points

29

#### More than simply a turbine

Through its combination of 20 % more power, less weight and the innovative Direct Stop Safety Technology, KaVo's MASTERtorque is a real alternative to high-speed instruments

Special colour editions anthracite and brown







## KaVo. Dental Excellence.